astron CLEARsplint®

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advantages of ASTRON CLEARsplint® that speak for themselves!

- extremely break resistant, hence optimal protection of the teeth
- thermoplastic flexibility which provides greater, comfortable and tension-free wearing comfort
- hypoallergenic, free of MMA (residual monomer content of MMA=0)
- self-adjusting with memory effect. This superb accuracy permits the closest fit
- always transparent - will not yellow with age, since amine-free

astron CLEARsplint®

the best solution for indications as easing the mandibular joint (bruxism), for the protection of implants, the adjustment of the occlusion etc.

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Attention storage of powder and liquid: store liquid in refrigerator at 5-7°C, powder at room temperature. But prior to mixing both, powder and liquid, must be refrigerated.
Please note: in order to avoid uncontrollable chemical reactions **please use only mixing bowls and -spatulas which do not get in contact with other resins.** Processing repairs as usual (refer to step 5 until 10 presspacking instructions). Prior to processing repairs, put some CLEARsplint® liquid on area to be repaired and let it dry for some seconds. Mixing ratio and temperature have to be respected precisely otherwise resin will not cure correctly.

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1. prepare and measure model and mark tooth equator - splint can be made above tooth equator
2. wax up splint with a wax plate - **Put the required amount of powder in the refrigerator (prior to mixing both components must be refrigerated)**
3. prepare protection wall and let it cure. Remove it and prepare on the tubers, resp. on both sides of the molar triangle an overflow canal
4. boil out the model and water it
5. isolate the still warm plaster model
6. **remove liquid from the refrigerator just prior to mixing and pouring**
7. mixing ratio: one part liquid to two parts powder
Important: respect mixing ratio thoroughly! For one splint approx. 8 - 10 g powder are required
8. mix it thoroughly **for approx. 5 - 10 seconds as quickly as possible and bubble-free** (use separate mixing cup and -spatula)
9. pour cold mixture immediately into the protection wall, reset carefully plaster model in the protection wall. Remove excess material
10. immerse in pressure pot as soon as possible at 2 bar and **50°C** for 45 minutes (**max. 53°C**)
11. remove denture mold from pressure pot and place it **under cold running tap water for 15 minutes** before removing splint from pour mold
12. trim and polish as usual
13. **immerse splint into handwarm water prior to placing it on the model or in the mouth (flexibilisation)**

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1. prepare and measure model and mark tooth equator - splint can be made above tooth equator
2. wax up splint with a wax plate and invest model into flask - **Put the required amount of powder in the refrigerator (prior to mixing both components must be refrigerated)**
3. prepare counter part and after curing heat flask in water bath
4. remove wax, boil out model and counter part, water and isolate them
5. **remove liquid from the refrigerator just prior to mixing**
6. mixing ratio: one part liquid to two parts powder
Important: respect mixing ratio thoroughly! For one splint approx. 8 - 10 g powder are required
7. mix thoroughly **for approx. 5 - 10 seconds as quickly as possible and bubble-free** (use separate mixing cup and spatula)
8. as soon as mix reaches a doughy state (5-8 minutes) pack into regular denture mold using normal press-packing procedure
9. immerse in **50°C** warm water and cure for 60 minutes at 2 bar (**max. 53°C**)
10. cool for **15 minutes under cold running tap water** before removing from mold
11. trim and polish as usual
12. **immerse splint into handwarm water prior to placing it on the model or in the mouth (flexibilisation)**

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